


Want to enhance transduction? RetroNectin is the answer!

Buy 1 ➡ 5%, Buy2 ➡ 10%, Buy3 ➡ 15%, Buy4 ➡ 20%

- ✓ Enhance viral transduction in hard to transfect cell types
- ✓ Works with Lenti-and retroviruses
- ✓ Less toxicity than other commonly enhancers like polybrene
- ✓ Increases expansion and naïve T-cell fraction of PBMCs
- ✓ Available also in GMP grade



Cat#	Product	Kit Size
T100B	RetroNectin®	2.5 mg
T110A	RetroNectin® Precoated Dish	10 Dishes
T100A	RetroNectin®	0.5 mg
T201	RetroNectin® GMP grade	2.5 mg

Please Contact Us

Danyel Biotech office:

08-9366066

mbc@danyel.co.il

Valid until: 31.12.17

